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© 2022 DH Tech. Dev. Pte. Ltd. RUO-MKT-02-14659-ZH-A

Determination of pyrethroids and macrocyclic lactone insecticides in spices and tea

Using the SCIEX 7500 system

Cathy Lane,¹ Sara Cheikh Ibrahim,² Tino Schroeder,² Susanne Hergett,² Jack Steed,¹ Jianru Stahl-Zeng,³ Roy Sperling²

¹SCIEX, UK; ²Bilacoon, Germany; ³SCIEX, Germany

Pyrethroids and macrocyclic lactones are groups of commonly used insecticides in the agriculture and horticulture industries. Macrocyclic lactones are naturally occurring, or semisynthetic, compounds produced as fermentation products in soil-dwelling *Streptomyces avermitilis*.¹ Pyrethroids, on the other hand, are synthetic, and were designed based on the naturally occurring family of pyrethrins, which were originally derived from chrysanthemum flowers.²

Due to the widespread use of these compounds in the environment, a comprehensive quantitative method is necessary to monitor and control their concentration in final food products destined for human consumption.

Here, a method has been developed using the SCIEX 7500 system for the simultaneous identification and quantification of pyrethroid and macrocyclic lactone insecticides at detection levels below the maximum residue level defined by the European Commission under regulation 2018/1514.³

Key features of the SCIEX 7500 system for the quantification of pesticides in spices and green tea

- Highly sensitive detection and quantification of avermectin (containing 96% avermectin B1a and 4% avermectin B1b), bifenthrin, cyfluthrin, cypermethrin, deltamethrin, fenvalerat, λ -cyhalothrin, milbemectin A3 and A4 and permethrin from QuEChERS extracts of spices and green tea
- Improved sensitivity over previous assays with lower limits of quantification (LLOQs) down to 0.02 ng/mL in solvent (Figure 1)
- Optimization of new parameter Q0D for milbemectin A3 in spices resulted in greatly reduced background and increased signal to noise
- Increased sensitivity allows for the use of lower sample injection volumes, increasing assay robustness and significantly reducing ion suppression in matrix


Figure 1. Signals obtained for avermectin B1a, bifenthrin and deltamethrin at their respective LLOQs. A good signal is obtained at a concentration level of 0.02 ng/mL, highlighting the sensitivity that can be achieved when using the SCIEX 7500 system.


Figure 2. Additional declustering of ions using Q0 dissociation. On the SCIEX 7500 system, additional declustering can be applied during LC-MS analysis using 2 modes: Q0D simple (left) and Q0D enhanced (right). Use of a voltage differential, applied either between the QJet ion guide rods and IQ0 lens or between the IQ0 lens and Q0 on the SCIEX 7500 system, can aid in breaking up clusters, removing interferences and increasing signal to noise.

Methods

Sample preparation: Avermectin (containing 96% avermectin B1a and 4% avermectin B1b), bifenthrin, cyfluthrin, cypermethrin, deltamethrin, fenvalerat, λ-cyhalothrin, milbemectin A3 and A4 and permethrin were diluted into mobile phase A for analysis. Calibration curves were constructed in mobile phase A by spiking in compounds at concentrations of 0.02–20 ng/mL.

Organic green tea (1 g) was milled to a fine powder, homogenized and added to 10 mL of water and 10 mL of acetonitrile. For the spices (organic paprika powder), a starting mass of 2 g was used. The sample was shaken for 1 min and added to Macherey-Nagel QuEChERS Mix I (ref 730970). The sample was shaken for another 3 min and centrifuged for 5 min. Next, 8 mL of the organic (upper) phase was added to Macherey-Nagel QuEChERS Mix III (ref 730648). The sample was shaken for 3 min, centrifuged for 5 min and 7 mL was then taken and acidified with 5% formic acid in acetonitrile (10 µL per 1 mL of sample taken). Compounds were spiked into the tea extract at 1 ng/mL and 10 ng/mL, corresponding to 10 ppb and 100 ppb in the original tea and 5 ppb and 50 ppb in the original spice sample. Samples were analyzed on the SCIEX 7500 system without further dilution.

Chromatography: Chromatographic separation was performed using the ExionLC AD system and a Phenomenex Synergi Fusion-RP (4 µm, 50 x 2.1 mm) column at a flow rate of 0.25 mL/min. Mobile phase A was 5 mM ammonium formate in 20% methanol and 80% water (v/v); mobile phase B was 5 mM ammonium formate in 90% methanol and 10% water (v/v). A gradient of 0–100% mobile phase B over 8 min was used. The total analysis time was 17 min, using either a 2 µL or 5 µL injection volume. The column temperature was 40°C.

Mass spectrometry: Multiple reaction monitoring (MRM) analysis was performed using a SCIEX 7500 system. The system was operated in positive electrospray ionization (ESI) mode using the OptiFlow Pro ion source. Data was acquired using SCIEX OS software. Source conditions, which were optimized using the guided optimization tool in SCIEX OS software, were ISV 3500 V, GS1 40, GS2 70, TEM 250°C and CUR 50.

Data processing: Data was processed using SCIEX OS software.


Figure 3. Sensitivity increase resulting from optimization of parameters for Q0D enhanced for milbemectin A3 in spices. XIC for milbemectin A3 using default Q0D value of -10 V (left); XIC for milbemectin A3 using optimized Q0D enhanced mode value of 25 V (right). Milbemectin was analyzed at a concentration of 50 ppb (2 µL injection volume). A signal to noise increase of 3.5-fold was observed after optimization of Q0D.

Method development and optimization

An MRM method was first developed and optimized on the SCIEX Triple Quad 5500 system and targeted more than 500 pesticides and their metabolites. Eleven of the most challenging compounds from the multi-component method were selected for evaluation on the SCIEX 7500 system, and the method was adapted. Using the guided optimization tool in SCIEX OS software, collision energies were re-optimized for the SCIEX 7500 system. Source conditions were also optimized using the guided optimization tool in SCIEX OS software.

Additional declustering of ions can be achieved on the SCIEX 7500 system by using a voltage differential applied either between the QJet ion guide and IQ0 lens (Q0 dissociation, Q0D simple) or between the IQ0 lens and Q0 rods (Q0D enhanced), illustrated in Figure 2. Q0D is useful for breaking up clusters, decreasing interferences and increasing signal to noise in some cases. Q0D was optimized for a subset of the analytes that were the most difficult to detect in green tea and spices. MRM extracted ion chromatograms (XICs) for milbemectin A3 in spices that were acquired using both default settings and optimized Q0D enhanced settings are shown in Figure 3. In this case, a signal to noise increase of 3.5-fold was observed after optimization of Q0D.

The final MRM assay comprised 1 MRM transition per compound. Separation of target compounds is shown in Figure 4.


Figure 4. Chromatographic separation of 10 compounds. Avermectin (containing 96% avermectin B1a and 4% avermectin B1b), bifenthrin, cyfluthrin, cypermethrin, deltamethrin, fenvalerat, λ -cyhalothrin, milbemectin A3 and A4 and permethrin at a concentration of 10 ng/mL (5 μ L injection volume) in mobile phase A.

Dilution series in solvent

The optimized MRM assay was applied to the analysis of a dilution series of target compounds in mobile phase A. Standards were prepared in the range 0.01–20 ng/mL; 5 μ L sample was injected per analysis. Each sample was analyzed 3 times. Examples of calibration curves and MRM XICs from the low end of the dilution series for 3 compounds are shown in Figure 5. LLOQs and associated statistics for all compounds are shown in Table 1.

Table 1. LLOQs and associated statistics for a dilution series of 10 compounds in solvent. Avermectin contains 96% avermectin B1a and 4% avermectin B1b. For each compound, linearity was observed from the LLOQ to the 20 ng/mL top standard. A linear regression with 1/x weighting was applied to the data.

| Compound | Q1 (m/z) | Q3 (m/z) | Retention time (min) | LLOQ (ng/mL) | # Values at LLOQ | % Accuracy at LLOQ | % CV at LLOQ | Linearity (r value) |
|------------------------|----------|----------|----------------------|--------------|------------------|--------------------|--------------|---------------------|
| AvermectinB1a | 890.5 | 305.1 | 9.35 | 0.02 | 3 of 3 | 142 | 7.4 | 0.9964 |
| AvermectinB1b | 876.5 | 291.1 | 9.12 | 0.50 | 3 of 3 | 99.9 | 6.6 | 0.9958 |
| Bifenthrin | 440.2 | 181.1 | 9.84 | 0.02 | 3 of 3 | 154 | 4.1 | 0.9937 |
| Cyfluthrin | 450.9 | 434.0 | 8.99 | 0.50 | 3 of 3 | 96.9 | 5.1 | 0.9890 |
| Cypermethrin | 433.1 | 191.0 | 9.09 | 0.02 | 3 of 3 | 106 | 10.5 | 0.9918 |
| Deltamethrin | 522.9 | 280.7 | 9.13 | 0.02 | 3 of 3 | 97.4 | 3.4 | 0.9924 |
| Fenvalerat | 437.0 | 125.0 | 9.22 | 0.20 | 3 of 3 | 93.4 | 5.6 | 0.9915 |
| λ -cyhalothrin | 467.0 | 225.0 | 9.06 | 0.05 | 3 of 3 | 110 | 12.6 | 0.9897 |
| Permethrin | 408.2 | 355.1 | 9.57 | 0.05 | 3 of 3 | 107 | 12.6 | 0.9914 |
| Milbemectin A3 | 546.1 | 511.2 | 9.24 | 0.05 | 3 of 3 | 99.3 | 6.6 | 0.9987 |
| Milbemectin A4 | 560.0 | 525.2 | 9.52 | 0.05 | 3 of 3 | 99.5 | 10.4 | 0.9965 |


Figure 5. MRM XICs from the low end of the calibration curves for avermectin B1a, bifenthrin and deltamethrin. An LLOQ of 0.02 ng/mL was obtained for each of these 3 compounds in solvent with linearity maintained across the entire dilution series (all r values ≥ 0.9924).

Compounds spiked into spices and green tea

Compounds were next spiked into QuEChERS extracts of spices and green tea at 1 ng/mL and 10 ng/mL, corresponding to 5 ppb and 50 ppb in spices and 10 ppb and 100 ppb in tea.

Quantification of target compounds in matrix against the dilution series in solvent was performed using both 2 μ L and 5 μ L injection volumes to assess the effect of lowering the injection volume on matrix ion suppression. It was found that the reduction in ion suppression observed for a 2 μ L injection compared with a 5 μ L injection mostly compensated for the reduction in injection volume to the extent that the sensitivity of the assay was not significantly impacted (Table 2 and Figure 6). Figure 7 demonstrates the background signal and ion suppression observed with a 2 μ L injection volume; however, the level of ion suppression was much less than the ion suppression observed with higher volume injections (Figure 6).

It is noted that for bifenthrin, an unexpected increase in signal was observed in the QuEChERS spiked samples compared with the dilution series in solvent. The source of this anomalous result is unknown and requires further investigation. See Table 3 and Table 4 for a summary of each compound's ion suppression results and concentrations of detection (mean calculated concentration) in the matrices analyzed.

| Compound | Q1 | Q3 | Spices, 2 μ L Injection | | | Tea, 2 μ L Injection | | |
|--------------------|-------|-------|-----------------------------|-------|--------|--------------------------|--------|---------|
| | | | Blank | 5 ppb | 50 ppb | Blank | 10 ppb | 100 ppb |
| AvermectinB1a | 890.5 | 305.1 | x | ✓ | ✓ | x | ✓ | ✓ |
| AvermectinB1b | 876.5 | 291.1 | x | x | x | x | x | x |
| Bifenthrin | 440.2 | 181.1 | ✓ | ✓ | ✓ | x | ✓ | ✓ |
| Cyfluthrin | 450.9 | 434.0 | x | x | ✓ | x | x | ✓ |
| Cypermethrin | 433.1 | 191.0 | x | x | ✓ | x | ✓ | ✓ |
| Deltamethrin | 522.9 | 280.7 | x | ✓ | ✓ | x | ✓ | ✓ |
| Fenvalerat | 437.0 | 125.0 | x | x | ✓ | x | ✓ | ✓ |
| Cyhalothrin,lambda | 467.0 | 225.0 | x | x | ✓ | x | ✓ | ✓ |
| Permethrin | 408.2 | 355.1 | x | x | ✓ | x | ✓ | ✓ |
| Milbemectin A3 | 546.1 | 511.2 | x | x | ✓ | x | x | ✓ |
| Milbemectin A4 | 560.0 | 525.2 | x | x | ✓ | x | x | ✓ |

| Compound | Q1 | Q3 | Spices, 5 μ L Injection | | | Tea, 5 μ L Injection | | |
|--------------------|-------|-------|-----------------------------|-------|--------|--------------------------|--------|---------|
| | | | Blank | 5 ppb | 50 ppb | Blank | 10 ppb | 100 ppb |
| AvermectinB1a | 890.5 | 305.1 | x | ✓ | ✓ | x | ✓ | ✓ |
| AvermectinB1b | 876.5 | 291.1 | x | x | x | x | x | x |
| Bifenthrin | 440.2 | 181.1 | ✓ | ✓ | ✓ | x | ✓ | ✓ |
| Cyfluthrin | 450.9 | 434.0 | x | x | ✓ | x | x | ✓ |
| Cypermethrin | 433.1 | 191.0 | x | x | ✓ | x | ✓ | ✓ |
| Deltamethrin | 522.9 | 280.7 | x | ✓ | ✓ | x | ✓ | ✓ |
| Fenvalerat | 437.0 | 125.0 | x | x | ✓ | x | ✓ | ✓ |
| Cyhalothrin,lambda | 467.0 | 225.0 | x | x | ✓ | x | ✓ | ✓ |
| Permethrin | 408.2 | 355.1 | x | x | ✓ | x | ✓ | ✓ |
| Milbemectin A3 | 546.1 | 511.2 | x | x | ✓ | x | x | ✓ |
| Milbemectin A4 | 560.0 | 525.2 | x | x | ✓ | x | x | ✓ |

Table 2. Detection of spices and green tea. Top: 2 μ L sample injected per analysis. Bottom: 5 μ L sample injected per analysis. No difference in sensitivity was observed when using a 2 μ L injection compared to a 5 μ L injection volume. This may be due to the reduction in ion suppression observed for a lower volume injection (see Figure 6).


Figure 6. Ion suppression in green tea at 100 ppb (top) and spices at 50 ppb (bottom) versus solvent for injection volumes of 2 μ L and 5 μ L. Each sample was analyzed 3 times at each injection volume, and concentrations were calculated against a dilution in solvent. The % ion suppression values were calculated based on expected concentrations. See Table 3 and Table 4 for more detail.


Figure 7. Analysis of 10 compounds in different matrices. MRM XICs for 10 compounds spiked into solvent (left), and into QuEChERS extracts of green tea (middle) and spices (right) using a 2 μ L injection volume. Data is from blank samples, shown at the top, and samples with a concentration of 10 ng/mL (100 ppb in green tea and 50 ppb in spices), shown at the bottom. High levels of background signals are present in the green tea and spice QuEChERS extracts compared with analysis of compounds spiked into solvent alone.

Table 3. Quantification of 10 compounds in green tea against a dilution series in solvent. A QuEChERS extract of green tea was spiked at 1 ng/mL and 10 ng/mL (10 ppb and 100 ppb) with each of 10 compounds and analyzed using the SCIEX 7500 system. Two injection volumes (2 μ L and 5 μ L) were used. Each sample was analyzed 3 times at each injection volume. Calculated concentrations against solvent curves were adjusted for injection volume.

| Compound | Solvent RT (min) | Green tea, 10 ppb, 2 μ L injection | | | Green tea, 100 ppb, 2 μ L injection | | | Green tea, 10 ppb, 5 μ L injection | | | Green tea, 100 ppb, 5 μ L injection | | |
|---|------------------|--|------------------------|--------------------------|---|------------------------|--------------------------|--|------------------------|--------------------------|---|------------------------|--------------------------|
| | | Matrix RT (min) | Mean calc. conc. (ppb) | Matrix % ion suppression | Matrix RT (min) | Mean calc. conc. (ppb) | Matrix % ion suppression | Matrix RT (min) | Mean calc. conc. (ppb) | Matrix % ion suppression | Matrix RT (min) | Mean calc. conc. (ppb) | Matrix % ion suppression |
| <i>AvermectinB1a</i> | 9.35 | 9.35 | 4.1 | 59.5 | 9.36 | 36.4 | 63.3 | 9.35 | 2.7 | 73.4 | 9.33 | 32.7 | 67.3 |
| <i>AvermectinB1b</i> | 9.12 | ND | ND | ND | ND | ND | ND | ND | ND | ND | ND | ND | ND |
| <i>Bifenthrin</i> | 9.84 | 9.84 | 23.9 | -138.8 | 9.85 | 264.0 | -163.9 | 9.84 | 14.3 | -42.5 | 9.83 | 164.0 | -63.8 |
| <i>Cyfluthrin</i> | 8.99 | ND | ND | ND | 9.02 | 39.8 | 60.3 | ND | ND | ND | 9.00 | 25.4 | 74.6 |
| <i>Cypermethrin</i> | 9.09 | 9.10 | 3.5 | 64.8 | 9.10 | 63.0 | 37.1 | 9.10 | 1.9 | 80.7 | 9.10 | 24.9 | 75.1 |
| <i>Deltamethrin</i> | 9.13 | 9.13 | 4.7 | 53.3 | 9.13 | 77.9 | 22.1 | 9.13 | 3.0 | 70.2 | 9.13 | 29.5 | 70.5 |
| <i>Fenvalerat</i> | 9.22 | 9.21 | 6.9 | 31.5 | 9.22 | 68.1 | 31.9 | 9.21 | 3.8 | 62.2 | 9.21 | 44.2 | 55.8 |
| <i>λ-cyhalothrin</i> | 9.06 | 9.05 | 3.5 | 65.5 | 9.07 | 56.3 | 43.7 | 9.05 | 2.5 | 74.9 | 9.05 | 24.6 | 75.4 |
| <i>Permethrin</i> | 9.57 | 9.56 | 7.4 | 26.0 | 9.57 | 96.5 | 3.5 | 9.56 | 9.8 | 1.9 | 9.56 | 72.1 | 27.9 |
| <i>Milbemectin A3</i> | 9.24 | ND | ND | ND | 9.23 | 47.2 | 52.9 | ND | ND | ND | 9.21 | 27.5 | 72.5 |
| <i>Milbemectin A4</i> | 9.52 | ND | ND | ND | 9.51 | 64.3 | 35.8 | ND | ND | ND | 9.45 | 32.2 | 67.9 |

ND = not detected

Table 4. Quantification of 10 compounds in spices against a dilution series in solvent. A QuEChERS extract of spices was spiked at 1 ng/mL and 10 ng/mL (5 ppb and 50 ppb) with each of 10 compounds and analyzed using the SCIEX 7500 system. Two injection volumes (2 μ L and 5 μ L) were used. Each sample was analyzed 3 times at each injection volume. Calculated concentrations against solvent curves were adjusted for injection volume.

| Compound | Solvent RT (min) | Spices, 5 ppb, 2 μ L injection | | | Spices, 50 ppb, 2 μ L injection | | | Spices, 5 ppb, 5 μ L injection | | | Spices, 50 ppb, 5 μ L injection | | |
|---|------------------|------------------------------------|------------------------|--------------------------|-------------------------------------|------------------------|--------------------------|------------------------------------|------------------------|--------------------------|-------------------------------------|------------------------|--------------------------|
| | | Matrix RT (min) | Mean calc. conc. (ppb) | Matrix % ion suppression | Matrix RT (min) | Mean calc. conc. (ppb) | Matrix % ion suppression | Matrix RT (min) | Mean calc. conc. (ppb) | Matrix % ion suppression | Matrix RT (min) | Mean calc. conc. (ppb) | Matrix % ion suppression |
| <i>AvermectinB1a</i> | 9.35 | 9.33 | 2.6 | 47.8 | 9.34 | 32.6 | 34.8 | 9.29 | 0.7 | 86.8 | 9.29 | 5.8 | 88.4 |
| <i>AvermectinB1b</i> | 9.12 | ND | ND | ND | ND | ND | ND | ND | ND | ND | ND | ND | ND |
| <i>Bifenthrin</i> | 9.84 | 9.82 | 8.8 | -75.3 | 9.84 | 69.0 | -38.0 | 9.81 | 1.5 | 70.6 | 9.82 | 15.4 | 69.2 |
| <i>Cyfluthrin</i> | 8.99 | ND | ND | ND | 9.02 | 20.0 | 60.2 | ND | ND | ND | 9.00 | 9.6 | 80.8 |
| <i>Cypermethrin</i> | 9.09 | ND | ND | ND | 9.09 | 23.6 | 52.9 | ND | ND | ND | 9.07 | 21.0 | 58.1 |
| <i>Deltamethrin</i> | 9.13 | 9.13 | 1.7 | 66.8 | 9.13 | 21.6 | 56.9 | 9.13 | 0.6 | 88.3 | 9.12 | 5.6 | 88.9 |
| <i>Fenvalerat</i> | 9.22 | ND | ND | ND | 9.21 | 28.9 | 42.3 | ND | ND | ND | 9.19 | 9.7 | 80.6 |
| <i>λ-cyhalothrin</i> | 9.06 | ND | ND | ND | 9.04 | 41.6 | 16.8 | ND | ND | ND | 9.02 | 13.3 | 73.4 |
| <i>Permethrin</i> | 9.57 | ND | ND | ND | 9.56 | 5.8 | 88.4 | ND | ND | ND | 9.54 | 1.4 | 97.2 |
| <i>Milbemectin A3</i> | 9.24 | ND | ND | ND | 9.21 | 22.3 | 55.5 | ND | ND | ND | 9.17 | 5.4 | 89.2 |
| <i>Milbemectin A4</i> | 9.52 | ND | ND | ND | 9.48 | 30.6 | 38.8 | ND | ND | ND | 9.49 | 6.7 | 86.7 |

ND = not detected

Conclusions

To summarize, the SCIEX 7500 system provides impressive levels of sensitivity, robustness and accuracy for trace quantification of insecticides in food matrices. In this study, excellent sensitivity has been demonstrated with LLOQ values down to 0.02 ng/mL. Quantification of pyrethroids and macrocyclic lactones in both green tea and spices has also been shown. By using a higher sensitivity LC-MS platform, injection volumes could be reduced, lowering ion suppression without impacting the overall assay sensitivity. The reduced amount of matrix injected for large studies can help to improve total system uptime.

References

1. ScienceDirect, [Avermectin](#).
2. Ccancapa-Cartagena, A.; Masiá, A.; Picó, Y. (2017) Simultaneous determination of pyrethroids and pyrethrins by dispersive liquid-liquid microextraction and liquid chromatography triple quadrupole mass spectrometry in environmental samples. [Anal Bioanal Chem 409\(20\), 4787–4799 \(2017\)](#).
3. European Commission, [Commission Regulation \(EU\) 2018/1514](#), October 2018.
4. Enabling new levels of quantification. [SCIEX technical note, RUO-MKT-02-11886-A](#).
5. Sensitivity gains for the evolution of routine bioanalysis. [SCIEX technical note, RUO-MKT-02-11885-A](#).
6. Tomlin, C. The pesticide manual – a world compendium, 13th ed.; [British Crop Protection Council \(Alton, Hampshire, UK\)](#), 2003.
7. European Commission, [Food Safety: EU Pesticides Database](#).
8. European Commission, [Food Safety: Maximum Residue Levels](#).

The SCIEX clinical diagnostic portfolio is For In Vitro Diagnostic Use. Rx Only. Product(s) not available in all countries. For information on availability, please contact your local sales representative or refer to www.sciex.com/diagnostics. All other products are For Research Use Only. Not for use in Diagnostic Procedures.

Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd. or their respective owners in the United States and/or certain other countries (see www.sciex.com/trademarks).

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Determination of 10 kinds of Neonicotinoid Pesticide Residues in Foods of Plant Origin by High Performance Liquid Chromatography-Mass Spectrometry

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- [1] 丁曉波, 張華, 劉世堯, 等. 柑橘果品營養學研究現狀[J]. 園藝學報
- [2] Syed Bilal Hussain, Cai-Yun Shi, et al., Recent Advances in the Regulation of Citric Acid Metabolism in Citrus Fruit. Critical Reviews in Plant Sciences
- [3] Baldwin, E. A. Citrus fruit. In: Biochemistry of Fruit Ripening (pp. 107–137).
- [4] 羅曾玲, 徐素珍. 柑橘的甜酸風味及生產管理措施分析; 羅曾玲 徐素珍. 鄉村科技
- [5] Qualitative flexibility combined with quantitative power -Using the SCIEX ZenoTOF 7600 LC-MS/MS system, powered by SCIEX OS software. SCIEX technical note, RUO-MKT-02-13053-A.

- [6] Etienne, A., Genard, M., et al., What controls fleshy fruit acidity? A review of malate and citrate accumulation in fruit cells. J. Exp. Bot. 64: 1451–1469.
- [7] 曾祥國. 不同種類和產區柑橘糖酸含量及組成研究〔D〕. 華中農業大學, 2005.
- [8] LIN Qiong, et al. Low temperature induced changes in citrate metabolism in Ponkan (Citrus reticulata Blanco cv · Ponkan) fruit during maturation · Plos One · 2016 ·

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© 2022 DH Tech. Dev. Pte. Ltd. RUO-MKT-02-14784-ZH-A

Quantification of Whey Protein Content in Infant Formula by CE-SDS

Gao Tie,¹ Chen Hongxu,¹ Wu Jianhua,¹ Feng Ping²

¹ CE & Biopharma, SCIEX, China

² Wyeth Nutrition, China

Introduction

Milk is almost always the first source of protein for human growth and development. Mammalian milk includes two major groups of proteins: whey proteins, a family of globular proteins, and casein proteins, a family of phosphoproteins. They are found in a roughly 60:40 ratio in human milk and in other ratios in other mammalian milks. Whey proteins are also a standard component of most infant formulas. According to the People's Republic of China National Food Standard GB 10765, milk-based infant formula refers to liquid or powdered products made mainly from milk and milk protein products; with addition of adequate amount of vitamins, minerals, and/or other ingredients; and produced and processed by physical methods only, that is suitable for normal infants.

Figure 2. PA 800 Plus Pharmaceutical Analysis System

To verify the claims of infant formula vendors and ensure the quality, safety, and efficacy of infant formula, it is essential to have an accurate method for measuring whey content. Previously, however, there was no reliable method to separate the two major protein types and quantify the whey protein


Figure 1. Electropherogram of infant formula.

content. The challenge was compounded by pasteurization, during which the Maillard reaction can cause proteins to form complexes with sugars. These complexes make effective separation of the various proteins even more difficult.

This note explores an effective, practical method based on capillary electrophoresis—sodium dodecyl sulfate (CE-SDS)—for separating and quantifying the whey protein content in infant formula. Analysis using the CE-based SCIEX PA 800 Series Pharmaceutical Analysis System clearly separated whey, casein proteins, and immunoglobulins (Figure 1), and provided accurate, reliable quantitation of the whey proteins. The method was validated across multiple laboratories and demonstrated to be reliable and reproducible.¹ This method was approved by the AOAC Expert Review Panel for Stakeholder Panel on Infant Formula and Adult Nutritionals (SPIFAN).¹

Key Features

- CE-SDS on the PA 800 Plus system provides excellent separation resolution, clearly separating the whey proteins, casein proteins, and immunoglobulins in both protein standard and infant formula samples
- Separations are relatively fast; under 30 minutes in the case of this analysis
- Validation across multiple laboratories found the method to be reliable and reproducible

Experimental

Reagents and Supplies

The SDS-MW Analysis Kit (SCIEX, P/N 390953) supplied the capillary, SDS gel buffer, 10 kDa marker, MW marker, and sample buffer. Contents are listed in Table 1. Other reagents and supplies included: β -mercaptoethanol (SIGMA P/N M7154 or M6250); milk protein standards (SIGMA: α -lactalbumin, L6010; β -lactoglobulin, L3908; IgG, I5506; BSA, A1933; β -casein, C6905; α S-casein, C6780; κ -casein, C0406. ALAR Ingredients: CGMP). Infant formula was provided by a local dairy manufacturer. Deionized water (Millipore).

| | Time (min) | Event | Value | Duration | Inlet vial | Outlet vial | Summary | Comments |
|---|------------|--------------------|----------|-----------|------------|-------------|---------------------------------------|---|
| 1 | | Rinse - Pressure | 20.0 psi | 10.00 min | BI:D1 | BO:D1 | forward | 0.1 N NaOH rinse to clean capillary surface |
| 2 | | Rinse - Pressure | 20.0 psi | 5.00 min | BI:E1 | BO:E1 | forward | 0.1 N HCl rinse to neutralize capillary surface silanol group |
| 3 | | Rinse - Pressure | 20.0 psi | 2.00 min | BI:F1 | BO:F1 | forward | ddH ₂ O rinse to remove the acid residue |
| 4 | | Rinse - Pressure | 70.0 psi | 10.00 min | BI:B1 | BO:B1 | forward | SDS Gel rise to fill the capillary |
| 5 | 0.00 | Separate - Voltage | 15.0 KV | 10.00 min | BI:C1 | BO:C1 | 5.00 Min ramp, reverse polarity, both | SDS Gel for voltage equilibration |
| 6 | | | | | | | | |

Figure 3. Column conditioning method.

Sample Preparation

Milk protein standard mixture

Milk proteins were weighed and dissolved in water according to the infant formula ratio.

Infant formula

Infant formula (100 mg \pm 4 mg) was added to 1.5 mL centrifuge tubes. Water (1 mL) was added to each tube. Tubes were vortexed and oscillated until the infant formula was fully dissolved.

Sample running pre-solution

Sample pre-running solution was prepared by mixing sample buffer and the 10 kDa Marker (internal standard), both from the SDS-MW Analysis Kit, in an 84:1 ratio.

Denaturation

10 μ L of each sample solution was pipetted into separate 1.5 mL microcentrifuge tubes. 85 μ L of sample running pre-solution and 5 μ L of β -mercaptoethanol were added to each microcentrifuge tube. Each tube was mixed well and then heated in a water bath at 100° \pm 5° C for 10 minutes. The tubes were allowed to cool to room temperature and then vortexed. Finally, each sample was transferred to a separate injection vial.

| Items | Part Number |
|--------------------------------|-------------|
| 50 μ m ID capillary column | 338451 |
| SDS Gel Buffer | A30341 |
| 10 kDa marker | A26487 |
| MW marker | A22196 |
| SDS sample buffer | |
| 0.1N NaOH | |
| 0.1N HCl | |

Table 1. Items in SDS-MW Analysis Kit

Separation and Detection

All separations were carried out on a capillary electrophoresis-based SCIEX PA 800 Plus Pharmaceutical Analysis System equipped with a photodiode array (PDA) detector operated at 220 nm, 2 Hz. Aperture was 2 (100 x 200 μm). The system was equipped with a fused-silica capillary column (50 μm ID, 30 cm total / 20 cm effective length). Sample storage and capillary temperatures were both 25° C. Conditioning and separation methods are show in Figures 3 and 4, respectively.

Data Analysis

Data analysis was performed using SCIEX 32 Karat Software.

Integration

Integration was carried out with Caesar integration off.

The electropherogram of the protein standard mixture was integrated. Migration of the whey and casein proteins are defined by the comparison to the E-grams of every single protein standard shown in Figure 5.

For the infant formula electropherogram, the autointegration start point was after the solvent peak and approximately 0.4 minutes before the 10 kDa marker peak. The auto integration end point was the end of the casein protein peaks group, the position of which was defined according to the standard mixture. Manual integration was performed for high MW whey proteins that eluted after the casein proteins.

| | Time (min) | Event | Value | Duration | Inlet vial | Outlet vial | Summary | Comments |
|----|------------|--------------------|----------|-----------|------------|-------------|--|--|
| 1 | | Rinse - Pressure | 70.0 psi | 3.00 min | BI:D1 | BO:D1 | forward, In / Out vial inc 8 | 0.1 N NaOH rinse to clean capillary surface - Automatic increment every 8 runs |
| 2 | | Rinse - Pressure | 70.0 psi | 1.00 min | BI:E1 | BO:E1 | forward, In / Out vial inc 8 | 0.1 N HCl rinse to neutralize capillary surface silanol group - Automatic increment every 8 runs |
| 3 | | Rinse - Pressure | 70.0 psi | 1.00 min | BI:F1 | BO:F1 | forward, In / Out vial inc 8 | Water rinse to remove the acid residue - Automatic increment every 8 runs |
| 4 | | Rinse - Pressure | 70.0 psi | 10.00 min | BI:B1 | BO:B1 | forward, In / Out vial inc 8 | SDS Gel rinse to fill the capillary with SDS gel - Automatic increment every 8 runs |
| 5 | | Wait | | 0.00 min | BI:A1 | BO:A1 | In / Out vial inc 8 | ddH2O, use for dipping to clean capillary tip - Automatic increment every 8 runs |
| 6 | | Wait | | 0.00 min | BI:A4 | BO:A4 | In / Out vial inc 8 | ddH2O, use for dipping to clean capillary tip - Automatic increment every 8 runs |
| 7 | | Inject - Voltage | 5.0 KV | 20.0 sec | SI:A1 | BO:C1 | Override, reverse polarity | Sample injection |
| 8 | | Wait | | 0.00 min | BI:B4 | BO:B4 | In / Out vial inc 8 | ddH2O, use for dipping to avoid sample carry over - Automatic increment every 8 runs |
| 9 | 0.00 | Separate - Voltage | 15.0 KV | 30.00 min | BI:C1 | BO:C1 | 1.00 Min ramp, reverse polarity, both, In / Out vial inc 8 | SDS Gel for separation - Automatic increment every 8 runs |
| 10 | 5.00 | Autozero | | | | | | |
| 11 | | | | | | | | |

Figure 4. Separation method.

Calculation

Whey protein content in the infant formula was calculated using the following equations:

$$Aw.c. = 1.4 \times Aw \quad (1)$$

where: Aw.c. = corrected integrated area of the whey components

Aw = total integrated area of the whey components

1.4 = correction factor

$$\text{Whey protein \%} = Aw.c. / (Aw.c. + Acn) \times 100 \quad (2)$$

where: Acn = integrated area of the casein components

The correction factor is used to adjust the response factor difference of whey proteins and casein proteins when detected at UV 220nm, since caseins contain more aromatic amino acids (t) which contribute higher response for the amount to area than whey proteins do.

Results and Discussion

Milk Protein Standard Mixture

The milk protein standard composed of whey proteins (α -lactalbumin, β -lactoglobulin, CGMP, IgG, and BSA) and casein proteins (β -casein, α S1-casein, α S2-casein, κ -casein) were analyzed by CE-SDS. All the casein proteins migrated between the IgG light chain and IgG heavy chain. The whey proteins and casein proteins were completely separated from each other (Figure 5).

Infant Formula

The electropherogram of the infant formula sample was previously shown in Figure 1. Whey protein content was calculated to be 61.97%, in line with the whey protein content of natural human milk.


Figure 5. Electropherogram of protein standard mix. The first integrated peak is the 10 kDa marker.

Conclusion

Capillary electrophoresis with SDS was successfully applied to the quantification of whey protein content in infant formula. The CE-SDS method fully separated the whey proteins, casein proteins, and immunoglobulins, facilitating quantification of the whey proteins. Validation in multiple laboratories has shown the method to be an accurate, reliable solution to the challenging problem of verifying whey content.

References

1. Ping Feng, Christophe Fuerer, Adrienne McMahon. Quantification of Whey Protein Content in infant Formulas by Sodium Dodecyl Sulfate-Capillary Gel Electrophoresis (SDS-CGE): Single-Laboratory Validation, First Action 2016.15; *Journal of AOAC International*, Vol. 100, No. 2, 2017

Who is SCIEX? SCIEX company's global leadership and world-class service and support in the capillary electrophoresis and liquid chromatography-mass spectrometry industry have made it a trusted partner to thousands of the scientists and lab analysts worldwide who are focused on basic research, drug discovery and development, food and environmental testing, forensics and clinical research.

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Identification and Quantitation of Sterols in rice by SCIEX Triple Quad™ 3500 system

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Determination of Irganox compounds extracted from food packaging using 4 food simulants

Quantification of Irganox chemicals using the QTRAP 6500+ system

Lakshmanan Deenadayalan¹, Sashank Pillai¹, Jessica Smith², Jack Steed² and Jianru Stahl-Zheng³
¹SCIEX, India; ²SCIEX, UK; ³SCIEX, Germany

Introduction

This technical note demonstrates the accurate and precise quantification of 16 Irganox compounds spiked into food simulants and in a plastic food packaging extract. Most Irganox compounds met the acceptable accuracy (70-130%) and precision criteria (%CV <15%) when quantifying from a single solvent external calibration curve, eliminating the need for matrix-matched standards. The lowest calibration standard measured ranged from 0.025 to 0.50 ng/mL, demonstrating the sensitivity of the QTRAP 6500+ system.

Human health can be impacted by chemicals that migrate into foods and beverages from food contact materials (FCMs), such as wrappers and containers. Understanding the risk of FCM chemicals requires sensitive and accurate analytical methods. Extractables are chemicals that are released when the FCMs are stressed, such as using solvents, elevated temperatures, solvent exposure time and agitation. Conversely, leachables are chemicals that migrate under ambient conditions.

The method presented in this technical note focuses on extractables found in food packaging. The method was developed for the quantification of a diverse suite of 16 Irganox compounds found in food packaging using 4 food simulants, as outlined by the EU regulations.¹

Key features of the QTRAP 6500+ system for the analysis of Irganox compounds

- Two external calibration curve preparations that were dependent on compound-specific solubility were used to analyze food packaging samples in 4 different food simulants, based on EU regulations
- The sensitivity of the QTRAP 6500+ system enabled accurate quantification below the 10 ng/mL EU regulatory level
- Linearity was observed between 0.025 and 50 ng/mL with a r^2 value >0.99
- Post-spike sample accuracy was observed between 70% and 130% with %RSD <15% when quantified using an external calibration curve (n=6)


Figure 1. Overlaid extracted ion chromatograms (XICs) of 16 Irganox compounds. Positive ion (top panel) and negative ion (bottom panel) analytes were acquired in electron spray ionization (ESI) mode and are shown for a 50 ng/mL standard.

Methods

Standard preparation: All standards were weighed to approximately 1 mg, dissolved in 1 mL of methanol and sonicated for 5 min. Irganox 1529 (n-hexane) and Irganox E201 (ethanol) were available as neat standard solutions. A 50 µg/mL mixed standard was prepared from the individual stock solutions with a final solvent composition of 10:90 (v/v), water/methanol.

Calibration curve preparation: Stock solutions were prepared in methanol and spiked into calibration solution 1 or 2, depending on compound-specific solubility. Calibration solutions were composed of 50:50 (v/v), water/ethanol with 3% acetic acid by volume (calibration curve 1) and 50:50 (v/v), water/ethanol (calibration curve 2). Thirteen of the 16 Irganox compounds were prepared as described for calibration curve 1. The final calibration curves ranged from 0.025 to 50 ng/mL.

Food simulant spikes: Compounds were spiked into each of the 4 food simulants to yield final in-sample concentrations of 0.10, 1.00 and 2.5 ng/mL but were then diluted using the appropriate solvent in a 1:2 ratio to ensure complete solubility (Table 1). Therefore, the final in-vial concentrations were 0.033, 0.33 and 0.83 ng/mL. The food simulant spike samples were quantified using both calibration curves 1 and 2.

Food packaging material extract preparation: A 1 g sample of a plastic food packaging material (yogurt container) was weighed into each 10 mL glass centrifuge tube and 10 mL of each simulant was added (Table 1). The tubes were vortexed for 5 minutes and incubated at 40°C for 24 hours. Samples were filtered through a 0.22 µm syringe filter and diluted using the appropriate dilution solvent in a 1:2 ratio (Table 1).

Post-extraction spike of food packaging material: Extracts were prepared as described above. The final extract was spiked with the mixed Irganox standard to yield an in-sample concentration of 10 ng/mL. The resulting solution was then diluted with the appropriate solvent in a 1:2 ratio (Table 1).

Table 1. Food simulants and dilution solvents used to evaluate plastic materials or articles intended for food contact based on the EU commission regulation (10/2011).

| Food simulant | Dilution solvent |
|----------------------|--|
| 3% Acetic acid (w/v) | 100% Ethanol |
| 10% Ethanol (w/v) | 50 % Ethanol with 3% acetic acid |
| 20% Ethanol (w/v) | 50 % Ethanol with 3% acetic acid |
| 50% Ethanol (w/v) | Samples were injected without dilution |

Chromatography: The ExionLC AD system was used with a Phenomenex Luna Omega Polar C18 analytical column (C18 100Å, 3 µm, 100 x 4.6 mm). In addition, an Agela Ghost Hunter column (50 x 4.6 mm) was used as a delay column to separate the analytical peak from the LC pump contamination. The LC gradient conditions used are shown in Table 2. Mobile phase A was water with 10mM ammonium formate and 0.1% formic acid by volume. Mobile phase B was methanol. The flow rate was 0.8 mL/min and the column oven temperature was 40°C.

Table 2. LC gradient used for the separation of 16 Irganox compounds.

| Time (min) | %A | %B |
|------------|----|----|
| 0.01 | 10 | 90 |
| 1.0 | 10 | 90 |
| 2.0 | 2 | 98 |
| 10.5 | 2 | 98 |
| 10.6 | 10 | 90 |
| 12.0 | 10 | 90 |

Mass spectrometry: Samples were analyzed using the SCIEX QTRAP 6500+ system with electrospray ionization (ESI) and polarity switching. Data were acquired using the Scheduled MRM algorithm. The target scan time was 0.80 sec for positive polarity and 0.30 sec for negative polarity. The optimized source and gas parameters and compound-specific MRM parameters used are shown in Tables 3 and 4, respectively.

Table 3. Optimized source and gas parameters.

| Parameter | Value |
|-------------------|--------------|
| Curtain Gas | 40 psi |
| CAD Gas | Medium |
| Ion Spray Voltage | 4500/-4500 V |
| Temperature | 500°C |
| GS1 | 90 psi |
| GS2 | 35 psi |

Data processing: All data were processed using SCIEX OS software, version 2.1.6. For the consistency with the EU regulations, concentrations are reported on an in-sample basis which is representative the method sensitivity, whereas the in-vial concentrations represent the instrument sensitivity.

Table 4. Optimized compound-specific parameters used for analysis.

| Compound | Q1 (m/z) | Q3 (m/z) | DP (V) | CE (V) | CXP (V) |
|-----------------|----------|----------|--------|--------|---------|
| Irganox PS 800 | 515.4 | 329.2 | 150 | 20 | 13 |
| Irganox 1010 | 1194.8 | 785.4 | 130 | 73 | 13 |
| Irganox 1330 | 792.6 | 219.1 | 125 | 37 | 12 |
| Irganox MD1024 | 553.4 | 237.1 | 120 | 38 | 10 |
| Irganox 565 | 589.3 | 250.0 | 160 | 62 | 15 |
| Irganox 5057 | 394.3 | 134.1 | 195 | 50 | 10 |
| Irganox 3125 | 1059.6 | 762.4 | 100 | 44 | 14 |
| Irganox 1098 | 637.5 | 581.4 | 130 | 32 | 11 |
| Irganox 259 | 656.6 | 415.3 | 80 | 34 | 10 |
| Irganox 1081 | 359.0 | 139.1 | 40 | 28 | 8 |
| Irganox E201 | 431.1 | 165.1 | 120 | 37 | 10 |
| Antioxidant 425 | 386.3 | 135.1 | 20 | 26 | 16 |
| Irganox 3114 | 801.5 | 784.5 | 88 | 20 | 12 |
| Irganox 1425 | 327.1 | 299.1 | -115 | -33 | -11 |
| Irganox 1035 | 641.3 | 423.2 | -130 | -42 | -13 |
| Irganox 1520 | 423.2 | 145 | -60 | -21 | -11 |

Chromatographic separation

Good separation was achieved for the 16 Irganox compounds using an optimized chromatography method (Figure 1). A fast polarity switching method was used to acquire data for both positive and negative ion compounds in a single injection. A delay column was introduced after the mixing chamber to reduce the background interference from the mobile phases, thereby improving the method detection limits.

Calibration for Irganox 1425_01: $y = 8.58881e4 x + 2...4$ ($r = 0.99927$, $r^2 = 0.99854$) (weighting: $1/x^2$)


Figure 2. Illustrative calibration curve for Irganox 1425. The calibration curve was prepared using 50:50 (v/v), water/ethanol with 3% acetic acid (calibration curve 1) as diluent over the range of 0.025-50 ng/mL. An r^2 value >0.99 was achieved.

Calibration curve sensitivity and linear dynamic range

Calibration curves were generated for all compounds analyzed across a concentration range of 0.025 to 50 ng/mL. Standards were prepared in the calibration 1 diluent (50:50 (v/v), water/ethanol with 3% acetic acid by volume), unless specified. Table 5 shows the accurate quantification achieved across this range with an r^2 value >0.99 . The lower limit of quantification (LLOQ) is represented by the lowest standard analyzed in the calibration range and this value varied from 0.025 to 0.5 ng/mL across Irganox compounds. Figure 2 shows an example calibration curve using the quantifier transition for Irganox 1425, which demonstrates the linear dynamic range of 0.025 to 50 ng/mL and $r^2 > 0.99$.

Table 5. Linearity range and correlation coefficient (r^2) achieved for the 16 Irganox compounds.

| Compound | Linearity range (ng/mL) | Correlation coefficient (r^2) |
|-----------------|-------------------------|-----------------------------------|
| Irganox PS 800 | 0.250–50 | 0.991 |
| Irganox 1010* | 0.500–50 | 0.993 |
| Irganox 1330 | 0.050–10 | 0.990 |
| Irganox MD1024 | 0.100–50 | 0.995 |
| Irganox 565* | 0.025–50 | 0.992 |
| Irganox 5057 | 0.050–25 | 0.993 |
| Irganox 3125 | 0.100–25 | 0.990 |
| Irganox 1098 | 0.100–50 | 0.999 |
| Irganox 259 | 0.025–50 | 0.996 |
| Irganox 1081 | 0.500–50 | 0.998 |
| Irganox E201 | 0.250–50 | 0.996 |
| Antioxidant 425 | 0.050–50 | 0.997 |
| Irganox 3114* | 0.500–50 | 0.992 |
| Irganox 1425 | 0.025–50 | 0.999 |
| Irganox 1035 | 0.025–50 | 0.991 |
| Irganox 1520 | 0.050–25 | 0.991 |

Note: All Irganox compounds were quantified using calibration curve 1 (50:50 (v/v), water/ethanol with 3% (v/v) acetic acid) except for those denoted with an asterisk, which were quantified using calibration curve 2 (50:50 (v/v), water/ethanol).

Accurate and precise quantification of Irganox compounds in food simulant spikes using a single solvent calibration curve

Excellent accuracy (70-130%) and precision (%CV <15%) were achieved for 13 of the 16 Irganox compounds in all 4 food simulant spikes when quantified with the calibration curve 1 solutions (Table 6). The remaining 3 compounds (Irganox 1010, 565 and 3114) showed good accuracy and precision with either calibration curve 1 or calibration curve 2 solutions, depending on the food simulant matrix. Excellent precision based on 6 replicates was achieved at the LLOQ of all compounds spiked in each food simulant. These results demonstrate that good data quality can be achieved for most Irganox compounds using a single solvent curve across the diversity of food simulant matrices specified in the EU regulations. Removing the need to run individual matrix-matched curves can save significant sample preparation and analysis time.

Food packaging material spikes

Overall, the 10 ng/mL post-extraction spikes of the representative plastic food material (for example, yogurt container) showed good accuracy among the 4 food simulants, indicating minimal quantitative bias due to the matrix (Figure 3). The individual data points in Figure 3 represent the average of 6 replicates per Irganox compound. The 10% and 20% ethanol food simulants showed 70-130% accuracy for all compounds except for Irganox 1010 in the 10% ethanol simulant. The 50% ethanol and 3% acetic acid extracts showed 3 and 4 Irganox compounds outside the 70-130% acceptable range, respectively. No Irganox compounds were detected in the yogurt container, prior to spiking.

It is recommended to analyze the Irganox compounds with accuracies that failed to meet the 70-130% accuracy criteria using matrix-matched calibration or with mass-labeled internal standards. Overall, these findings demonstrate that the majority of the Irganox compounds in the 4 food simulants can be accurately quantified without the need for internal standards.


Figure 3. Average accuracy (n=6) of post-extraction spikes from the plastic yogurt container. Final extracts were spiked to yield 10 ng/mL in-sample and then diluted 2:1 with the compatible diluent solvent. Dots represent averages for individual Irganox compounds and the curve width represents approximate frequency of data points.

Table 6. Average accuracy and %CV (n=6) for all 16 Irganox compounds at the lowest measured concentration. Concentrations are reported on an in-sample basis.

| Simulant | 3% Acetic acid | | 10 % Ethanol | | 20 % Ethanol | | 50% Ethanol | |
|------------------------|---------------------------------|---------------------|---------------------------------|---------------------|---------------------------------|---------------------|---------------------------------|---------------------|
| | In-sample concentration (ng/mL) | Avg. accuracy (%CV) | In-sample concentration (ng/mL) | Avg. accuracy (%CV) | In-sample concentration (ng/mL) | Avg. accuracy (%CV) | In-sample concentration (ng/mL) | Avg. accuracy (%CV) |
| Irganox PS 800 | 1.00 | 119 (11) | 1.00 | 120 (13) | 1.00 | 124 (13) | 1.00 | 94 (5.6) |
| Irganox 1330 | 1.00 | 114 (3.0) | 1.00 | 121 (2.8) | 1.00 | 120 (3.1) | 0.025 | 112 (3.3) |
| Irganox MD1024 | 0.10 | 113 (1.9) | 0.10 | 108 (6.8) | 0.10 | 109 (2.7) | 0.025 | 100 (5.9) |
| Irganox 5057 | 0.10 | 115 (3.7) | 0.10 | 117 (8.4) | 0.10 | 122 (4.0) | 0.025 | 122 (9.1) |
| Irganox 1098 | 0.10 | 118 (12) | 0.10 | 115 (7.9) | 0.10 | 107 (1.5) | 0.025 | 100 (1.8) |
| Irganox 259 | 0.10 | 114 (2.4) | 0.10 | 111 (2.2) | 0.10 | 110 (2.4) | 0.025 | 108 (0.8) |
| Irganox 1081 | 1.00 | 105 (13) | 1.00 | 101 (11) | 1.00 | 112 (11) | 1.00 | 105 (5.0) |
| Irganox E201 | 1.00 | 103 (4.8) | 2.50 | 113 (3.1) | 2.50 | 113 (9.7) | 1.00 | 125 (3.4) |
| Antioxidant 425 | 0.10 | 95 (13) | 0.10 | 88 (6.4) | 0.10 | 86 (12) | 0.025 | 70 (3.3) |
| Irganox 1425 | 1.00 | 108 (4.1) | 1.00 | 112 (2.9) | 1.00 | 108 (3.4) | 0.10 | 79 (11) |
| Irganox 1035 | 0.10 | 111 (9.6) | 0.10 | 109 (11) | 1.00* | 95 (5.3) | 0.025 | 106 (13) |
| Irganox 3125 | 1.00 | 106 (6.2) | 1.00 | 104 (5.6) | 1.00 | 103 (5.6) | 0.10 | 91 (6.9) |
| Irganox 1520 | 1.00 | 117 (3.0) | 1.00 | 103 (2.5) | 1.00 | 105 (2.3) | 0.025 | 117 (6.0) |
| Irganox 1010 | 2.50* | 106* (8.4) | 2.50 | 106 (12) | 2.50 | 96 (2.4) | 1.00* | 118* (7.3) |
| Irganox 565 | 1.00 | 96 (7.9) | 1.00* | 95* (2.6) | 1.00* | 85* (7.4) | 0.025* | 90* (5.9) |
| Irganox 3114 | 1.00* | 114* (5.4) | 1.00* | 121* (2.6) | 1.00* | 117* (3.3) | 1.00* | 89* (2.9) |

Note: All Irganox compounds were quantified using calibration curve 1 (50:50 (v/v), water/ethanol with 3% (v/v) acetic acid) except for those denoted with an asterisk, which were quantified using calibration curve 2 (50:50 (v/v), water/ethanol).

Conclusions

- A fast polarity switching method was developed to analyze 16 Irganox compounds in a single method
- The SCIEX QTRAP 6500+ system provided excellent sensitivity, reproducibility and linearity when quantifying Irganox compounds in different food simulants and post-extraction spikes in a plastic food container
- Two solvent-based external calibration curves were used to evaluate Irganox compounds in different simulants. Matrix-matched calibration curves were not required.
- R² values >0.99 were achieved for all compounds with %CV <15% across 6 replicates at the lowest measured concentration in all simulants, meeting the EU regulation requirements

References

1. Commission regulation (EU) No 10/2011 of 14th January 2011 on plastic materials and articles intended to come into contact with food ([EUR-Lex - 32011R0010 - EN - EUR-Lex \(europa.eu\)](#)).

The SCIEX clinical diagnostic portfolio is For In Vitro Diagnostic Use. Rx Only. Product(s) not available in all countries. For information on availability, please contact your local sales representative or refer to <https://sciex.com/diagnostics>. All other products are For Research Use Only. Not for use in Diagnostic Procedures.

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Application of the QTRAP® system in the screening of extractable substance of packaging materials

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Chen xi, Zhao xianglong, Li lijun, Guo lihai

SCIEX China

Key words: SCIEX OS; Packaging materials; Screening

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For Research Use Only. Not for use in Diagnostics Procedures.

AB Sciex is operating as SCIEX.

© 2019. AB Sciex. The trademarks mentioned herein are the property of AB Sciex Pte.

Ltd. or their respective owners. AB SCIEX™ is being used under license.

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A Rapid Quantitative Method of LC-MS /MS for 12 kinds of phenolic compounds in Food packaging material

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National Reference Laboratory for Food Contact Material (Changzhou)²

Keywords: SCIEX Triple Quad; Food packaging material; phenolic compound

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© 2022 DH Tech. Dev. Pte. Ltd. RUO-MKT-02-14599-ZH-A